

A summary of “Accurate Quantitation of Heroin Metabolite, 6-Monoacetylmorphine in Urine using IMCSzyme® for Hydrolysis”

Overview

Hydrolysis of drug metabolite conjugates is common practice to analyze drugs of abuse in human matrices on LC-MS/MS. β -glucuronidase enzyme is generally used to hydrolyze glucuronide conjugates. However, in the case of heroin use, where 6-monoacetylmorphine (6-MAM) may be present in urine, the conversion of 6-MAM to morphine may occur during hydrolysis. This study was performed to investigate the 6-MAM conversion by four commercially available β -glucuronidase enzymes.

Materials and Methods

All drug standards were purchased from Cerilliant Corporation. β -Glucuronidase enzymes were from Integrated Micro-Chromatography Systems, LLC (IMCSzyme®), Campbell Science (*Helix pomatia* and *Patella vulgata*). 200 μ L of urine was hydrolyzed with 100 μ L of enzyme solution (containing hydrolysis buffer and internal standard) at 55 °C with shaking. The hydrolysis times were 1 and 2 hours for each enzyme. The hydrolyzed sample was extracted with RP-S tip, provided by DPX Labs, LLC, and then analyzed by LC-MS/MS.

Results

In the urine sample hydrolyzed with IMCSzyme, there was less than 10 ng/mL, or less than 1% increase in morphine concentration and no significant difference in 6-MAM concentration over the 2 hour time incubation period. Comparing to other enzymes, the levels of morphine detected increased over the hydrolysis time, corresponding with the decrease in 6-MAM. The enzyme that converts 6-MAM into morphine the highest to the lowest is in the order: *Helix pomatia*, *Patella vulgata*, *Haliotis rufescens*, and IMCSzyme (Figure 1)

Conclusions

Due to the possible conversion of 6-MAM, most laboratories have to treat 6-MAM positive samples differently than other opiate positive samples. The results showed that while the three other commercially available enzymes significantly converts 6-MAM to morphine, IMCSzyme recombinant enzyme did not. The decrease in conversion of 6-MAM to morphine would allow for more accurate analytical runs of heroin positive samples even in the presence of a glucuronidase.

References

1. Bogusz, M., Maier, R., Driessen, S. Journal of Analytical Toxicology.1997; 21:346-355
2. Cone, E., Welch, P. Journal of Analytical Toxicology. 1991; 15: 1-7.
3. Smith, M., Shimomura, E., Summers, J., Paul, B. Journal of Analytical Toxicology. 2001; 25: 504-514
4. Boerner, U., Abbott, S., Roe, R. Drug Metabolism Reviews. 1975; 4(1): 39-73.
5. Mitchell, J., Paul, B. Journal of Analytical Toxicology.1991; 15(2): 49-53.
6. Yeh, S., Gorodetzky, C., McQuin, R. Journal of Pharmacology and Experimental Therapeutics.1976; 196(2): 249-256
7. Zuccaro, P., Ricciarello, R., Pichini, S., Pacifici, R., Altieri, I., Pellegrini, M. Ascenzo, G. Journal of Analytical Toxicology.1997; 21(4): 268-277
8. Cone, E., Jufer, R., Darwin, W. Journal of Analytical Toxicology. 1996; 20(6):379-392.
9. Aderjan, R., Skopp, G. Therapeutic Drug Monitoring. 1998; 20(5): 561-569. 10. ACD/ChemSketch, version 12, Advanced Chemistry Development, Inc., Toronto, On, Canada, www.acdlabs.com, 2014.




Figure 1. Concentration, in ng/mL, of 6-MAM and morphine after 0, 1, and 2 hours of incubation with the respective enzyme in previously drug free urine spiked with 6-MAM.