Using quantitative phosphoproteomics to understand key phosphorylation signaling pathways in HCT116 cells after chemotherapy drug treatment

<u>Todd Mullis^{1*}, Andrew Lee², Rebekah Woolsey³, David Quilici³, Qian Wang¹</u>

¹University of South Carolina, Department of Chemistry and Biochemistry, Columbia, SC; ²Integrated Micro-Chromatography Systems, Inc, Irmo, SC

³Mitch Hitchcock, Ph.D. Nevada Proteomics Center, 1664 N Virginia Street MS 0330, Reno, NV 89557

INTRODUCTION

We previously established an automated phosphopeptide enrichment method using an automated liquid handler and IMCStips[®]. Using this method, we

generated nearly 10,000 phosphopeptide IDs without prior fractionation from cell digest. We applied this automated method to probe the effect of a thymidylate synthase inhibitor, 5-fluorouracil (5-FU), on a common colorectal cancer (CRC) cell line, HCT116. 5-FU is a well-documented chemotherapeutic agent that is commonly used as a first-line treatment of many cancers, including CRC. However, innate and acquired resistance to the drug is common^{1,2}. There is a need for new drug targets that increase the drug's efficacy and prognostic markers that can predict the patient's sensitivity to the drug. In this study, comparative phosphoproteomic analysis using automated sample preparations is presented as an initial proof of concept to dissect the differentia phosphorylation events in HCT116 cells treated with 5-FU. The resulting data allows for the mapping of signaling cascades and monitoring of cellular events which may be useful in identifying new drug targets and predicting patient's sensitivity to the drug treatment.



MATERIALS AND METHODS

HCT116 cells were grown in Gibco™ DMEM/F-12 containing 15 mM of HEPES and 10% fetal bovine serum. Freshly thawed cells were passaged by trypsin/EDTA treatment at 80-90% confluence. Once the cells were passaged five times, they were treated with various concentrations of 5-FU for 48 hours. The treated cells were assayed to determine cell viability using the Promega™ CellTiter-Blue™ Cell Viability Assay. Stable isotope labelling with amino acids in cell culture (SILAC) was used to isotopically label HCT116 cells. They were grown in Thermo Scientific[™] DMEM:F-12 Media containing 10% dialyzed FBS. Mass spectrometry confirmed that incorporation of the heavy labeled amino acids (¹³C₆ L-Lysine and ¹³C₆¹⁵N₄ L-Arginine) was above 95%. Upon complete labelling, "light" cells were treated with 10 μN 5-FU and "heavy" cells were treated with 0.05% DMSO for 48 hours. The cells were lysed with RIPA buffer containing a protease and phosphatase inhibitor cocktail. Based on protein concentrations determined using the bicinchoninic acid (BCA) assay, the protein lysates were combined 1:1. The proteins were reduced in 10 mM TCEP at 56°C for 30 minutes and were then alkylated in 25 mM iodoacetamide at room temperature in the dark for 30 minutes. Following reduction and alkylation, the proteins were trypsin digested (1:50 enzyme to protein ratio) overnight at 37°C. For automated sample processing we used the Microlab® STAR™ from Hamilton Robotics. Five mg RP 300 μL IMCStips[®] were used for desalting prior to phosphopeptide enrichment with 10 mg PolyTi or 5 mg ZrO₂ 300 μL IMCStips[®].

- ...Orbitrap Fusion, Thermo Scientific Mass spectrometer
- .Ultimate 3000 RSLCnano, Thermo
- ...0.1% formic acid in acetonitrile • Mobile phase B.
- ..2% 90% B for 180 min LC gradient.
- ..40-cm in-house packed column with 1.8 μm Sepax GP-C18 (120Å, Sepax Technologies, Newark, DE).
- Column oven temperature...50°C

For global phosphopeptide identification, we used an Orbitrap Fusion mass spectrometer coupled to an UltiMate 3000 RSLCnano system. The mass spectrometer was operated in data-dependent mode at top speed. The dynamic exclusion duration was set to 60 s with a 10 ppm mass tolerance. Spectra were searched using the Sequest HT algorithm within Proteome Discoverer v2.2 (Thermo Scientific) in combination with the human UniProt protein FASTA database (2017 March, 20,119 entries). Search parameters were as follows: FT-trap instrument, parent mass error tolerance of 10 ppm, fragment mass error tolerance of 0.6 Da (monoisotopic), variable modifications of 15.995 Da (oxidation) on methionine and 79.966 Da (phosphorylation) on serine, threonine and tyrosine, fixed



Scheme 1. Experimental workflow of SILAC-based proteomic and phosphoproteomic analysis of 5-FU treated HCT116 cells. The workflow incorporates an automated desalting and phosphopeptide enrichment step that significantly reduces the hands-on time required for sample generation.

modification of 57.021 Da (carbamidomethylation) on cysteine. For SILAC quantitation, the heavy labelled cells were used as the control. The resulting data was analyzed using Ingenuity Pathway Analysis software. The overall scheme of the sample preparation and analysis process can be seen in **Scheme 1**.

RESULTS

We established an in-house 5-FU treatment of HCT116 cells by determining the 50% growth inhibition concentration (GI₅₀) for a 48-hour treatment. HCT116 cells were treated with various concentrations of 5-FU including, 1.56, 3.12, 6.25, 12.5, 25, 50, 100, 200, 400 and 800 µM for 48 hours. Using the CTB assay, the relative cell viability was determined for each 5-FU concentration (Table 1). A Monod equation was generated to calculate the GI 50 based on the data points from the CTB assay **(Figure 1)**. Based on this equation, the GI_{ro} for our HCT116 cell line was 42 μM. Previous studies have shown $1.064^{0.187}$ the GI₅₀ 5-FU concentration for the HCT116 cell line to be between 0.5-77 μ M^{3,4.} Using these results, 10 μ M 5-FU RFU_[5-FU] = RFUo _____ was selected as the concentration to be used for the cell treatment with a calculated theoretical growth inhibition of 35%. This established 5-FU treatment of HCT116 cells was then used to perform the SILAC-based quantitative

phosphoproteomic analysis.

| [5-FU], μM | 590 nm, RFU | Percent Inhibition, % |
|---------------|----------------|--------------------------|
| 0 | 3470.7 | 0% |
| 1.5625 | 3379.3 | 3% |
| 3.125 | 2689.8 | 22% |
| 6.25 | 2488.2 | 28% |
| 12.5 | 2012.6 | 42% |
| 25 | 1660.7 | 52% |
| 50 | 1586.9 | 54% |
| 100 | 1434.2 | 59% |
| 200 | 1469.0 | 58% |
| 400 | 1491.6 | 57% |
| 800 | 885.9 | 74% |

The relative growth inhibition percentage of HCT116 after 48-hour treatment with specified concentrations of 5-FU. The growth inhibition was calculated based the CTB assay. Percent inhibition was calculated using the control (0 µM 5-FU) as the relative value for 0% inhibition



Figure 1. Cell viability plot of HCT116 cells after 5-FU treatment at each concentration. The cell viability was quantified by 590 nm RFU values after incubation with CTB reagent. The displayed trendline is based on the generated Monod equation using the 590 nm RFU values from the CTB assay.



In previous studies, we established a fully automated desalting and phosphopeptide enrichment method using IMCStips[®] on a Hamilton Microlab[®] STAR[™], which allowed us to identify over 10,000 phosphopeptides from 200 µg of HCT116 cell lysate without fractionation. Here, as a proof of concept study, we used the established automated method to desalt and enrich for phosphopeptides from 200 μg of 1:1 mixed Heavy:Light cell lysate to quantitate protein fold changes in 5-FU treated HCT116 cells. For the analysis, 8 samples were desalted using RP IMCStips[®]. Two of the desalted samples were collected and vacuum dried without further enrichment. Six samples were further enriched using phosphopeptide enrichment IMCStips[®], 3 using PolyTi IMCStips[®] and 3 using ZrO₂ IMCStips[®]. Our previous studies showed a slight differential enrichment between the two resin types, so both PolyTi and ZrO₂ resin types were used for enrichment to obtain a comprehensive phosphorylation profile.

From the analysis, a total of 5,985 unique proteins and 67,250 unique peptides were identified from our analyzed samples. For the phosphopeptide enriched samples, the PolyTi IMCStips® samples averaged a total of 10,465 phosphopeptides identified and an average specificity of 85%. ZrO, IMCStips[®] samples averaged a total of 8,351 phosphopeptides identified with an average specificity of 90% (Figure 2). The unenriched samples processed using RP IMCStips[®] samples show very low phosphopeptide identifications, which highlights the need for the phosphopeptide enrichment.



Figure 2. A comparison of the total amou phosphopeptides and non-phosphopeptides ider by either PolyTi or ZrO, resin. All samples were processe with IMCStips[®] using our fully automated phosphopeptic enrichment method (n=3)

Proteins, #

2838

1862

enriched for phosphopeptides.

Peptides,

19785

18911

6788

6678

6422

4919

5175

4921

Table Further, the reproducibili[†] of the quantitation based or lds,# lds.# peptide fold change ratios were compared for both unenriched and enriched samples by plotting the log, fold change value for biochemical replicates 12388 (Figure 3). By Pearson PolvTi correlation, the r-values for the comparison o calculated fold change ratios for two replicates replicates from RP, PolyTi ZrO 2407 9138 and ZrO₂ were 0.88, 0.80

and 0.78 respectively. Based

on these high r-values, the sample preparation and analysis process is able to generate reproducible quantitative proteomic and phosphoproteomic data.



Reproducibility of Desalting Extraction- RP Peptides



Reproducibility of Phosphopeptid Enrichment - PolyTi Peptides

10365



PolyTi Sample 1- Log₂(Fold Change)

Reproducibility of Phosphopeptide Enrichment - ZrO, Peptides

Summary of generated MS data from 8 samples analyzed. 8 200 µg 1:1 Heavy:Light

enriched for phosphopeptides, whereas samples prepared with PolyTi or ZrO, were

HCT116 cell digest samples were analyzed. Samples prepared with RP resin were not

Reproducibility of Phosphopeptide Enrichment- ZrO₂ Peptides



Figure 3. A comparison peptide fold change ratios (log2) between two biochemical replicates of 200 µg of 1:1 mixed SILAC samples that were prepared using (A) the desalting extraction using RP IMCStips[®], (B) the desalting extraction using RP IMCStips[®] followed by the phosphopeptide enrichment using PolyTi IMCStips[®] or (C) the desalting extraction using RP IMCStips[®] followed by the phosphopeptide enrichment using ZrO₂ IMCStips[®]. R-values equal to 0.88, 0.80 and 0.78 respectively by Pearson correlation (P-Values < 0.00001).



Using the calculated fold change ratios, we then determined the number of proteins "significantly" decreased or increased for each sample analyzed. A protein was considered "significantly" decreased or increased if the calculated fold change was above 2 or below 0.5 with a calculated p-value of 0.05. For the unenriched samples an average of 685 ± 6 proteins were significantly decreased, and an average of 523 ± 18 proteins were significantly increased. For the phosphopeptide enriched samples, an average of 659 ± 20 proteins were significantly decreased and an average of 576 ± 27 proteins were significantly increased. Volcano plots for each sample were generated to show the overall distribution of protein fold changes.

Figure 4. Volcano plots of quantified proteins from each sample analyzed. The volcano plot is based on the proteins fold change and its associated p-value. (A) RP- Sample 1 Proteins: 679 Significantly Decreased, 5 Significantly Increased. (B) RP- Sample 2: 690 ntly Decreased, 541 Significantly Increas Decreased. 576 Significantly Increased. (D) P 31 Significantly Increased. (G) ZrO₂- Sample 2 Proteins 640 Significantly Decreased, 553 Significantly Increased. (H) ZrO,-Sample 3 Proteins: 652 Significantly Decreased, 584 Significantly Increased

To further determine the biological relevance of the data, Ingenuity Pathway Analysis (IPA) software was used. For this proof-of-concept study, we were primarily interested in analyzing the phosphopeptide enriched data set to understand the phosphorylation signaling pathways affected by the 5-FU treatment. The expression fold change data from the phosphopeptide enriched samples was analyzed using IPA to determine the affected signaling pathways (Figure 5).

From this analysis, the EIF2 signaling pathway was shown to be significantly inhibited. The EIF2 signaling pathway is generally associated with cellular growth and proliferation. It is known that in response to stress and amino acid starvation phosphorylation of eIF2a can inhibit translational initiation and by association cellular growth ^{1, 5, 6}. In our data set, we observed a collective decrease in kinases within the EIF2 pathway, which are generally associated with growth factors such as PI3K, AKT, c-RAF, potentially in response cell stress (Figure 6).

Further, a tox list was generated to understand the pathways associated with toxicity that were enriched in the analyzed dataset (Figure 7). The NRF2-Mediated Oxidative stress response pathway was shown to be increased within our dataset. This pathway has been implicated as a primary response pathway to 5-FU treatment due to an increase in



Figure 5. Canonical pathways identified using the proteins from the phosphopeptid enriched data set. The p-value represents the probability of random association of molecules with the identified pathway. A -log(P-value) threshold of 5 for this dataset. A positive z-score indicates an increase of the pathway, while a negative z-score indicates an inhibition of the pathway.

reactive oxygen species⁷. In previous studies, it has been shown that upregulation of NRF2 can lead to resistance to 5-FU treatment. Additionally, phosphorylation is thought to play a key role in the upregulation of NRF2 and NRF2-associated proteins⁸. The p53 signaling pathway was also shown to be affected in our dataset. It has recently been reported that p53 signaling pathway is crucial for initiating cell death in response to 5-FU treatment. In cells resistant to 5-FU treatment, the WNT-signaling pathway can inactivate the p53 signaling pathway cell death response².



© 2000-2010 QIACEN, All rights reserved

Figure 6. EIF2 signaling pathway with proteins identified and quantified within the phosphopeptide enriched sample set. The z-score for this pathway was calculated to be -2.46 indicating inhibition of the EIF2 signaling pathway. Molecules with purple bands indicate identification of the protein within the dataset.



Figure 7. Tox list generated by IPA software based on molecules from the phosphopeptide enriched dataset. The p-value with the identified pathway. A -log(P-value) threshold of 3 for this dataset.



The phosphopeptide enriched data set was further filtered by selecting proteins that were identified across all analyzed phosphopeptide samples, and were calculated to have a fold change of greater than 2 or less than 0.5. In total 22 proteins were identified, 11 proteins with a fold change greater than 2 and 11 proteins with a fold change less than 0.5 (Table 3). Lamin-B1 is of particular interest as it was recently reported that the overexpression of the protein leads to worse clinical outcomes in patients with colon cancer due to an increase in mitotic catastrophe allowing cancers cells to avoid apoptosis⁹. Additionally, upregulation of deoxycytidine kinase, an important enzyme in DNA synthesis, has been shown to correlate well with 5-FU sensitivity and has been highlighted as a potential prognostic marker for the drug treatment ¹⁰.

| Table 3 | | | |
|------------|-------------------------------------------------------------------|---------------------|-----------------|
| Uniprot Id | Protein Description | Average Fold Change | Phospho-Site(s) |
| Q8IV50 | LysM and putative peptidoglycan-binding domain-containg protein 2 | 100.00 | S24 |
| Q9Y285 | PhenylalaninetRNA ligase alpha subunit | 100.00 | S301 |
| Q8N6N3 | UPF0690 protein C1orf52 | 11.91 | T155 |
| Q92878 | DNA repair protein RAD50 | 9.77 | S635; T690 |
| Q92547 | DNA topoisomerase 2-binding protein 1 | 5.46 | S350; S888 |
| P27707 | Deoxycytidine kinase | 4.35 | S74 |
| Q8N4S0 | Coiled-coil domain-containing protein 82 | 4.03 | S282 |
| 014545 | TRAF-type zinc finger domain-containing protein 1 | 3.50 | S327; S415 |
| Q8WWA1 | Transmembrane protein 40 | 2.88 | S141 |
| Q08174 | Protocadherin-1 | 2.78 | S962; S984 |
| Q9Y6W5 | Wiskott-Aldrich syndrome protein family member 2 | 2.35 | S293 |
| Q5VTL8 | Pre-mRNA-splicing factor 38B | 0.44 | S529 |
| Q5JVS0 | Intracelluar hyaluronan-binding protein 4 | 0.36 | S108 |
| P22234 | Multifunctional protein 4 | 0.31 | S27 |
| 043583 | Density-regulated protein | 0.30 | S73 |
| Q6AI08 | HEAT repeat-containing protein 6 | 0.30 | S337; S643 |
| Q15717 | ELAV-like protein 1 | 0.29 | S202 |
| Q9UKX7 | Nuclear pore complex protein Nup50 | 0.28 | S296 |
| P20700 | Lamin-B1 | 0.26 | S23; S391; T575 |
| Q9H910 | Hematological and neurological expressed 1-like protein | 0.24 | S97 |
| Q9UN79 | Transcription factor SOX-13 | 0.01 | S334 |
| Q86U06 | Portable RNA-binding protein 23 | 0.01 | S149 |

Based on these results, this study shows the usefulness of SILAC-based phosphoproteomic studies for understanding drug treatments affects, and for identifying phosphorylated proteins that could be further studied to determine if they could be used as prognostic markers for drug sensitivity or resistance. In future studies we plan to expand upon our current dataset by performing a phospho-tyrosine enrichment in addition to the global phosphopeptide enrichment as our current dataset is primarily enriched for serine- and threonine-phosphorylated peptides.

CONCLUSION

In this proof of concept study, we used an established, fully automated desalting and phosphopeptide enrichment method to conduct a SILAC-based phosphoproteomic study of a 5-FU chemotherapy treatment of HCT116 cells. In this study we were able to identify 5,985 unique proteins and 67,250 unique peptides. From the mass spectrometry analyses, we quantified fold changes of phosphoproteins in response to the chemotherapy treatment and correlated those changes with changes in biological signaling pathways. Additionally, we identified 22 phosphoproteins that were significantly increased or decreased in all phosphopeptide enriched samples in response to the 5-FU treatment. These proteins could be explored further as potential drug targets or prognostic markers for 5-FU treatments.

REFERENCES

1. R. B. Hamanaka, B. S. Bennett, S. B. Cullinan and J. A. Diehl, *Mol Biol Cell*, 2005, **16**, 5493-5501.

- 2. L. He, H. Zhu, S. Zhou, T. Wu, H. Wu, H. Yang, H. Mao, C. SekharKathera, A. Janardhan, A. M. Edick, A. Zhang, Z. Hu, F. Pan and Z. Guo, Exp Mol Med, 2018, 50, 101. 3. K. Bracht, A. M. Nicholls, Y. Liu and W. F. Bodmer, *Br J Cancer*, 2010, **103**, 340-346.
- 4. P. M. De Angelis, D. H. Svendsrud, K. L. Kravik and T. Stokke, *Mol Cancer*, 2006, **5**, 20.
- 5. D. Zhou, L. R. Palam, L. Jiang, J. Narasimhan, K. A. Staschke and R. C. Wek, *J Biol Chem*, 2008, **283**, 7064-7073.
- 6. C. Kaehler, J. Isensee, T. Hucho, H. Lehrach and S. Krobitsch, *Nucleic Acids Res*, 2014, **42**, 6436-6447.
- 7. U. Ozer, K. W. Barbour, S. A. Clinton and F. G. Berger, *Mol Pharmacol*, 2015, **88**, 970-981.
- 8. A. Giudice, C. Arra and M. C. Turco, *Methods Mol Biol*, 2010, 647, 37-74
- 9. M. Izdebska, M. Gagat and A. Grzanka, *Int J Oncol*, 2018, **52**, 89-102.
- 10. F. McAllister, et al., *Cancer Biol Ther*, 2014, **15**, 688-698.

*Contact: Todd Mullis – mullist@email.sc.edu Artwork and layout created by IMCS, Inc. © 2019 IMCS, Inc. All rights reserved. IMCStip is a registered trademark of Integrated Micro Chromatography Systems, Inc.

Interested in more information about my poster?

Scan the QR code. I will follow up with you by email after the conference.

