



Automation of Sample Preparation and Buffer Exchange for Multi-Attribute Method



P. Nikki Sitasuwan¹; Thomas W. Powers²; Tiffany Medwid²; Yuko Ogata³; Nancy S. Nightlinger³; Richard S Rogers⁴; Casey Snodgrass⁵; Pamela M. Quizon¹; L. Andrew Lee¹
¹Integrated Micro-chromatography Systems, Inc, Irmo, SC; ²Pfizer Inc., Chesterfield, MO; ³Just - Evotec Biologics, Seattle, WA; ⁴Bristol-Myers Squibb, Seattle, Washington; ⁵Hamilton Company, Reno, NV




Figure 1. Sample preparation workflow using either manual or automated preparation. Experiments were performed at two independent sites using antibody stock with concentrations ranging from 1 – 10 mg/mL.

INTRODUCTION

- Multi-Attribute Method (MAM) – extension of peptide mapping applied to the characterization of a biopharmaceutical; leverages advances in ultrahigh-performance liquid chromatography coupled to mass spectrometry (UHPLC-MS) and data processing software.
- The combination of these technologies via MAM allows the simultaneous detection, identification, quantitation, and monitoring of molecular attributes of biopharmaceuticals.
- Current sample preparation for MAM relies on manual buffer exchange to remove excess denaturant and accelerate subsequent trypsin digestion efficiency.

GOALS

- To compare manual sample preparation using Bio-Spin P6 spin columns vs. automated sample preparation with SizeX IMCStips on a Hamilton Microlab STAR
- To evaluate robustness and reproducibility of the automated method by testing at two independent sites
- To collect precision data from three separate timepoints to determine consistency of the automated sample preparation method

MATERIALS AND METHODS

- SizeX IMCStips[®] were provided by IMCS. For manual buffer exchange, Bio-Spin[®] P6 columns (Bio-Rad) were used.
- Antibody stocks (1 to 10 mg/mL) were denatured, reduced, and alkylated
- The denatured protein samples (0.25 to 1 mg/mL) buffer exchanged using Bio-Spin P6 columns (manual) or SizeX IMCStips (automated)
- Desalted antibody was digested with trypsin to generate peptides for MAM analysis. Peptide samples analyzed on Thermo Q Exactive plus; data processed in BioPharma Finder. Known critical quality attributes were quantified and new peptide peaks were screened.

RESULTS

SITE ONE




Figure 2. Total ion chromatograms (TICs) of peptide samples from manual and automated methods using 1 mg/mL antibody stock.




Figure 3. Fc glycosylation patterns of automated samples compared to manual samples.




Figure 4. Modifications observed in automated samples compared to manual samples (Deam: deamidation, Oxi: oxidation).

SITE TWO




Figure 5. Antibody recovery from manual and automated MAM using Bio-Spin P6 columns and SizeX IMCStips, respectively (n = 4).




Figure 6. Total ion chromatograms (TICs) comparing automated (SizeX IMCStips) and manual preparation (Bio-Spin P6 Columns).




Figure 7. Comparison of modification percentages of samples prepared using Bio-Spin P6 columns (manual), and Hamilton STAR with SizeX IMCStips (automated). (A) Deamidation, succinimide and methionine oxidation. (B) Proline amidation and C-terminal lysine. (C) Fc glycosylation patterns (n = 4).




Figure 8. Precision of automated MAM sample preparation performed at three different timepoints with 4 replicates at each timepoint.

DISCUSSION & CONCLUSION

Automated MAM sample preparation utilizing SizeX IMCStips on Hamilton STAR liquid handling showed comparable precision and improved reproducibility over manual preparation. Automating tedious and repetitive sample preparation is a promising improvement for obtaining accurate and reproducible data for monitoring critical quality attributes of biopharmaceuticals. Furthermore, it could also serve as a platform to systematically optimize preparation conditions.

ACKNOWLEDGEMENTS

The authors would like to thank the MAM Consortium for their guidance.