

INTRODUCTION

- Urine drug testing is a common practice for monitoring compliance of prescribed opioid medications.
- Screening and confirmatory assays, such as liquid chromatography coupled with mass spectrometry (LC-MS/MS), benefit from the use of β -glucuronidases (β GUS) that deconjugate glucuronidated analytes.
- In this study we compared two recombinant β GUS, IMCS's IMCSzyme RT and Kura Biotech's B-One in two ways:
 - Ninety-six clinical samples were hydrolyzed with either IMCSzyme RT or B-One using 20 μ g of protein in the reaction.
 - Control and pooled clinical samples positive for morphine, oxymorphone, codeine and hydromorphone were hydrolyzed with multiple protein quantities ranging from 0 – 200 μ g of protein in the reaction.

MATERIALS AND METHODS

Drug free urine control (UTAK) was fortified with glucuronide standards equivalent to 100 ng/mL of free base. Opioid-positive urine specimens were obtained from a national testing laboratory. Control and pooled clinical samples were hydrolyzed with either IMCSzyme RT (IMCS) or B-One (Kura Biotech) with multiple protein concentrations ranging from 0 – 200 μ g (Table 1). Protein concentration for B-One was based on Bradford assay using BSA as standards, resulting in approximately 1 μ g/ μ L. IMCSzyme RT is quantified by A280, resulting in 2 μ g/ μ L. Hydrolyzed samples were diluted to 40% methanol with 5% formic acid in methanol elution solvent and filtered through a β -Gone plus plate to remove protein (Phenomenex). Samples were then diluted to 10% methanol with water prior to injection onto LC-MS/MS.

Table 1. IMCSzyme RT and B-One reaction set up.

IMCSzyme® RT					
Urine Sample (μ L)	IMCSzyme® RT at 2 mg/mL (μ L)	IMCSzyme® RT in reaction (μ g)	Water (μ L)*	Room Temperature Hydrolysis Buffer (μ L)	Internal Standard (μ L)
100	0	0	50	300	20
	2	4	48		
	4	8	46		
	6	12	44		
	8	16	42		
	10	20	40		
	30	60	20		
	50	100	0		

*water was added after hydrolysis so every sample had equivalent final volume

B-One®					
Urine Sample (μ L)	B-One® at 1 mg/mL (μ L)	B-One® in reaction (μ g)	Water (μ L)*	Hydrolysis Buffer (μ L)	Internal Standard (μ L)
100	0	0	200	B-One® comes all-in-one with hydrolysis buffer	20
	20	20	180		
	60	60	140		
	100	100	100		
	200	200	0		

*water was added after hydrolysis so every sample had equivalent final volume

10 μ L of diluted sample was injected on a Thermo Scientific™ Vanquish™ UHPLC system coupled with a Thermo Scientific™ Endura™ Triple Quadrupole Mass Spectrometer (Table 2). IMCSzyme RT and B-One had separate calibration curves using drug free urine as the matrix. All calibration curves had $r^2 \geq 0.992$ and calibration controls were within $\pm 20\%$. LLOQ is 10 ng/mL.

Table 2. HPLC and MS Parameters

Analysis Time	10 minutes
Column	Phenomenex Kinetex® 2.6 μ m Biphenyl 100 Å, 4.6 x 50 mm
Column Temperature	40°C
Mobile Phase A	0.1% formic acid in water
Mobile Phase B	0.1% formic acid in acetonitrile
Electrospray voltage	1000 V
Sheath Gas	55 Arb
Auxiliary gas	11 Arb
Sweep gas	1 Arb
Ion Transfer Tube Temperature	300°C
Vaporizer Temperature	300°C

RESULTS

IMCSzyme RT and B-One hydrolysis of morphine, oxymorphone, codeine and hydromorphone in control or clinical samples at multiple protein concentrations.




Figure 1. (A) Recovery of morphine across different amounts of protein in the sample indicates that 8 μ g of IMCSzyme RT was sufficient in control set to achieve > 80% recovery, in both the control \bullet and pooled clinical sample Δ . (B) The corresponding glucuronide peak diminishes as more protein is added, but < 20% glucuronide remains within the first sample of 4 μ g of protein for both control \bullet and pooled clinical sample Δ . (C) Similar recovery plot with B-One as increasing amounts of protein is added indicates more than IMCSzyme RT is needed to achieve similar recoveries in the control \bullet , and greater amount of B-One is needed in the pooled clinical sample Δ . (D) The corresponding glucuronide remaining in the sample supports this need for more protein in pooled clinical sample Δ in comparison to control sample \bullet . B-One protein to volume conversion is 1 μ g to 1 μ L, and IMCSzyme RT conversion is 2 μ g to 1 μ L.




Figure 3. Recovery of codeine and decreased detection of codeine-glucuronide across increasing amounts of protein. (A) Recovery of codeine with IMCSzyme RT for both control \bullet and pooled clinical sample Δ are closely aligned with over 80% recovery starting around 16 μ g of protein. (B) Corresponding codeine glucuronide diminishes with increase in IMCSzyme RT in comparable manner for control \bullet and clinical samples Δ . (C) Codeine recovery is significantly challenged with B-One, showing 200 μ g addition to achieve complete conversion in control sample \bullet but below 80% recovery threshold in clinical sample Δ . (D) Corresponding codeine glucuronide samples continue to be detected in both control \bullet and pooled clinical sample Δ , which suggests incomplete recovery of codeine. B-One protein to volume conversion is 1 μ g to 1 μ L, and IMCSzyme RT conversion is 2 μ g to 1 μ L.




Figure 2. Recovery of oxymorphone and corresponding decrease in glucuronidated metabolite across increasing amounts of protein. (A) Both control \bullet that was fortified with glucuronide standards equivalent of 100 ng/mL of free base and pooled clinical sample Δ performs in a comparable manner with consistent recoveries of oxymorphone, and (B) corresponding glucuronides diminish with increase amounts of protein in comparable manner for IMCSzyme RT in control \bullet and pooled clinical sample Δ . (C) Recovery of oxymorphone in control \bullet exhibits a classic plot whereas the pooled sample Δ shifts the recovery to the right, suggestive of sensitivity of the protein towards competitive inhibitors in the pooled sample, and (D) similar corresponding levels of glucuronide are detected in control \bullet and pooled clinical sample Δ . B-One protein to volume conversion is 1 μ g to 1 μ L, and IMCSzyme RT conversion is 2 μ g to 1 μ L.




Figure 4. Recovery of hydromorphone and corresponding glucuronide across increasing protein amounts. (A) Comparable recoveries are observed for control \bullet and pooled sample Δ , and (B) for glucuronide detection decreases as the corresponding free base increases in control \bullet and pooled clinical sample Δ . (C) While more protein is needed to achieve the similar recovery as IMCSzyme RT, this recoveries in control \bullet and pooled urine Δ are comparable, and (D) corresponding decrease in hydromorphone glucuronide decreases as greater amounts of protein is added to the sample. B-One protein to volume conversion is 1 μ g to 1 μ L, and IMCSzyme RT conversion is 2 μ g to 1 μ L.

Due to limited supply of B-One, ninety-six clinical samples were hydrolyzed with either IMCSzyme RT or B-One using 20 μ g of protein in the reaction containing 100 μ L of sample, 300 μ L of RTB for IMCSzyme RT or no buffer for B-One for fixed time of 15 minutes at room temperature (20-25 °C). All samples contained 20 μ L internal standard. Each sample set had its own calibration standards. Samples are considered positive (+) with ≥ 25 ng/mL of drug and negative (-) with < 25 ng/mL of drug. Results for Oxymorphone and Codeine are shown in Table 3 and 4.

		B-One Hydrolysis	
		(+)	(-)
IMCSzyme RT Hydrolysis	(+)	20	16
	(-)	0	60

Table 3. Oxymorphone

In Agreement:

- 20 samples as both positive
- 60 samples as both negative

In Disagreement:

- None of the samples that were positive by B-One were negative by IMCSzyme RT
- 16 samples disagreed where they were positive by IMCSzyme RT but negative by B-One

		B-One Hydrolysis	
		(+)	(-)
IMCSzyme RT Hydrolysis	(+)	0	5
	(-)	0	91

Table 4. Codeine

In Agreement:

- 91 samples as both negative

In Disagreement:

- None of the samples that were positive by B-One were negative by IMCSzyme RT
- 5 samples disagreed where they were positive by IMCSzyme RT but negative by B-One

CONCLUSION

- Clinical samples contain endogenous urine chemicals that can affect hydrolysis efficiencies likely due to presence of denaturing chemicals like urea, and other competing glucuronidated metabolites.
- In specified testing parameters, both DFU and pooled urine samples required more B-One than IMCSzyme RT to achieve > 80% recovery for all four opiates.
- When using 20 μ g of protein and 100 μ L of urine with 15 minute room temperature incubation, there were more samples below the 25 ng/mL threshold for samples processed with B-One than samples processed with IMCSzyme RT.
- Using less β GUS product may result in higher number of false negatives.
- Future studies may entail comparing larger amount of B-One or modified buffer ratios to ensure appropriate pH conditions are achieved.

REFERENCES

- Lee LA, McGee AC, Sitasuwan P, Tomashek JJ, Riley C, Muñoz-Muñoz AC, Andrade L. Factors Compromising Glucuronidase Performance in Urine Drug Testing Potentially Resulting in False Negatives. *J Anal Toxicol.* 2022 Jul 14;46(6):689-696. doi: 10.1093/jat/bkab090.
- Foster FD & Stuff JR. (2020). *Automated Hydrolysis Extraction and Analysis of Synthetic Cathinones in Urine using a Robotic Autosampler and LC-MS/MS Platform.* [Online]. Kura Biotech Finden B-One. <https://www.kurabiotech.com/finden-b-one>. Accessed March 30, 2023.
- Technical DataSheet B-One. [Online]. Kura Biotech Finden B-One. <https://www.kurabiotech.com/finden-b-one>. Accessed March 30, 2023.
- Test Catalog – CSMTU: Controlled Substance Monitoring Targeted Profile, 17 Drug Classes, Mass Spectrometry, Random, Urine. [Online]. Mayo Clinic Laboratories. <https://www.mayocliniclabs.com/test-catalog/overview/615292>. Accessed March 30, 2023.
- Rozas M, Callejas JL, Berner C, Vargas M, Wallach E, Andrews J, Novoa C, inventors; Kura Biotech Inc, assignee. One-step formulation for all-in-one hydrolysis solution. United States patent application US 17/324,067. 2021 Nov 25.

DOWNLOAD
POSTER



Scan this QR code
to get access to
this poster online.